



# Nuance® PowerScribe® 360 Critical Results

## Critical Results Desktop Installation Using MSI Package

1. Login to <https://ps360criticalresults.com> as an administrator and click on Tools – Downloads.




2. Click and download the **MSI PACKAGE FOR PUSH INSTALLATIONS**, which downloads the PS360CRDesktop.msi file.



3. If you have a previous version of CR Desktop application installed on the workstations, first run the following command from command prompt, to uninstall the application in quiet mode:

```
msiexec.exe /x "{E19A083B-BC30-4CE3-8F6E-617E0483463E}" /qb
```

# Nuance® PowerScribe® 360 Critical Results



```
Administrator: C:\Windows\system32\cmd.exe
C:\>msiserver /x "<E19A083B-BC30-4CE3-8F6E-617E0483463E>" /qb
C:\>
```

4. The generic command to install CR Desktop using MSI file:

```
msiexec.exe /passive /log "<LOGFILEPATH>" /i "<PS360CrDesktop MSI PATH>"  
SiteType=<SITETYPE>
```

where SiteType=1 for Production and SiteType=2 for Test (Staging)

5. In the example below, the MSI file is located under D:\Software and D: drive is where the Desktopinstall.log will be written. Change the path as required. You can also use a shared drive path for the MSI package, depending on the permissions to access it.

```
msiexec.exe /passive /log "D:\DesktopInstall.log" /i "D:\Software\PS360CRDesktop.msi" SiteType=1
```



```
Administrator: C:\Windows\system32\cmd.exe
C:\>msiexec.exe /passive /log "D:\DesktopInstall.log" /i "D:\Software\PS360CRDesktop.msi" SiteType=1
```

6. To use a Proxy Server, following parameters can be added:

```
PROXYIP  
PROXYPORT  
PROXYUSERNAME  
PROXYPASSWORD
```

Example:

```
msiexec.exe /passive /log "D:\DesktopInstall.log" /i "D:\Software\PS360CRDesktop.msi"  
SiteType=1 PROXYIP=192.168.1.1 PROXYPORT=5075 PROXYUSERNAME=admin  
PROXYPASSWORD=admin123
```

Nuance®  
**PowerScribe® 360**  
**Critical Results**

7. The application should be installed under C:\Program Files (x86)\Nuance\Critical Results Desktop on a Windows 64 bit workstation and C:\Program Files\Nuance\Critical Results Desktop on Windows 32 bit workstation.