

# Cerner RadNet Integration Quick Reference Card

## PowerScribe® 360 | Reporting (version 2.0 and higher) Optional Feature

**NOTE:** Nuance employees are not authorized to adjust or configure the settings on a customer's PACS. The customer, or the customer's PACS representative, must make any necessary changes to the PACS configuration or settings.

### Feature Overview

**NOTE:** This integration operates in **Slave Mode only**.

- Can launch a *PowerScribe 360 | Reporting* report using a PACS work list
- Can be installed on a different PC than the PACS/RIS
- Can be installed on the same PC as the PACS/RIS

#### Multiple Sites

If the system has multiple *sites* and the integration type is master, duplicate the master mode integration under each site. If the integration type is slave do not set up duplicate integrations under each site. Set up all slave integrations under the primary site (first site created). When a user logs in, no matter what site they are working under, the system will read the slave integrations created under the primary site and confirm the workstation is set up with the correct directory structures.

### Site PACS Settings

This API integration operates as a Slave and **requires no configuration**. The integration works through the .NET Remoting port (9090), labeled as the **IRadWhere port** in *PowerScribe 360 | Reporting*.

In addition, make sure that the **Broadcast IRadWhere events** check box is selected.

Both of these settings are located in the PowerScribe 360 | Reporting admin portal under **Setup > System > Configuration** on the **PowerScribe 360 Client** tab.

The screenshot shows the 'System Configuration' window for the 'PowerScribe 360 Client'. The 'PowerScribe 360 Client' tab is active. The following settings are visible:

- Dragon user directory: http://MLB-PS360-25/DragonUsers
- Dragon share username: psservice
- Dragon share password: \*\*\*\*\*
- Dragon maximum compute speed: 70%
- IRadWhere port: 9090** (highlighted with a red rectangle)
- Broadcast IRadWhere events: ☒** (highlighted with a red rectangle)
- Days to keep application logs: 10 days

Nuance® PowerScribe® 360

This page intentionally left blank.