

Intelrad IntelViewer Integration Quick Reference Card

PowerScribe® 360 | Reporting (version 2.0 and higher) Optional Feature

NOTE: Nuance employees are not authorized to adjust or configure the settings on a customer's PACS. The customer, or the customer's PACS representative, must make any necessary changes to the PACS configuration or settings.

Feature Overview

NOTE: Version Requirement: Your PACs should be at IV-4-6-1-R16 (4-6-1-P129) or -IV-4-7-1-R5 (4-7-1-P83) or higher.

In **Master** mode:

- Application launch
- Single login/logout
- Launch PACS images using *PowerScribe 360 | Reporting* work list
- Can be installed on the same PC as the PACS/RIS

In **Slave** mode:

- Mark exam as dictated/read
- Launch *PowerScribe 360 | Reporting* report using PACS work list
- Can be installed on the same PC as the PACS/RIS

Multiple Sites

If the system has multiple *sites* and the integration type is master, duplicate the master mode integration under each site. If the integration type is slave do not set up duplicate integrations under each site. Set up all slave integrations under the primary site (first site created). When a user logs in, no matter what site they are working under, the system will read the slave integrations created under the primary site and confirm the workstation is set up with the correct directory structures.

Site PACS Settings: Master

1. Go to **Setup > Sites > PACS** section.
2. Click the **Edit**  icon for the PACS you want to configure as Master (or click the **Add**  icon to create a new PACS entry).
3. From the **Type** drop-down list, select **Intelrad**.
4. In the **Name** field, type **Intelrad**. (The **Description** field is optional.)
5. In the **Mode** field, select **Master**.
6. In the **Master** text field, enter the following four parameters (**which must be provided by Intelrad or the customer**):
 - Intelrad server URL
 - Intelviewer password
 - Intelviewer fallback user
 - Authentication tokenExample: **https://integration.intelrad.com;nuance;nuancefallback;0d855d6732104911c9a6712d90a2cc3c**
7. Click the **Save**  icon. Your entry should look like the example illustration below.

PACS:	Type	Name	Description	Mode
		Intelrad	Intelrad	Master:https://integration.intelrad.com;nuance;nuancefallback;0d855d6732104911c9a6712d90a2cc3c
				

Page 1

Site PACS Settings: Slave

When using *PowerScribe 360 | Reporting* in the **Slave** mode, **you do not need to add an entry in the PACS settings section.**

This page intentionally left blank.