

Carestream (Kodak) Integration Quick Reference Card

PowerScribe® 360 | Reporting (version 2.0 and higher) Optional Feature

NOTE: Nuance employees are not authorized to adjust or configure the settings on a customer's PACS. The customer, or the customer's PACS representative, must make any necessary changes to the PACS configuration or settings.

Site PACS Settings

This API integration operates as a Slave or Master with no configuration required. The integration works through the .NET Remoting port (9090).

Multiple Sites

If the system has multiple *sites* and the integration type is master, duplicate the master mode integration under each site. If the integration type is slave do not set up duplicate integrations under each site. Set up all slave integrations under the primary site (first site created). When a user logs in, no matter what site they are working under, the system will read the slave integrations created under the primary site and confirm the workstation is set up with the correct directory structures.

This page intentionally left blank.