

Capitalization Tab

Use this tab to set capitalization rules for the names of persons, departments, drugs, and so on.

Nuance PowerScribe® 360 | Reporting: Speech Formatting - Windows Internet Explorer

Speech Formatting - Last Update: 9/12/2013 9:34:03 AM [Reset to Defaults](#) | [Save and Close](#)
| [Close Window](#)

General **Capitalization** Numbers, Units, and Dates Abbreviations Section Headings Verbs Miscellaneous

Person name capitalization:	Standard	[Standard]
Department name capitalization:	Standard	[Standard]
Drug name capitalization:	Standard	[Standard]
Capitalize the word after a new line:	<input checked="" type="checkbox"/>	[True]
Capitalize "allergy" and "allergies" as "ALLERGY" and "ALLERGIES":	<input type="checkbox"/>	[False]
Capitalize "against medical advice" as "AGAINST MEDICAL ADVICE":	<input checked="" type="checkbox"/>	[True]
Expand and capitalize "AMA" as "AGAINST MEDICAL ADVICE":	<input type="checkbox"/>	[False]
Capitalize "code blue" as "CODE BLUE":	<input checked="" type="checkbox"/>	[True]
Capitalize "code R" as "CODE R":	<input checked="" type="checkbox"/>	[True]
Capitalize "code status" as "CODE STATUS":	<input checked="" type="checkbox"/>	[True]
Capitalize "code S" as "CODE S":	<input checked="" type="checkbox"/>	[True]
Capitalize "code T" as "CODE T":	<input checked="" type="checkbox"/>	[True]
Capitalize "do not intubate" and "do not resuscitate" as "DO NOT INTUBATE" and "DO NOT RESUSCITATE":	<input checked="" type="checkbox"/>	[True]
Expand and capitalize "DNI" and "DNR" as "DO NOT INTUBATE" and "DO NOT RESUSCITATE":	<input type="checkbox"/>	[False]
Capitalize "durable power of attorney" as "DURABLE POWER OF ATTORNEY":	<input checked="" type="checkbox"/>	[True]
Capitalize "fish" as "FISH":	<input type="checkbox"/>	[False]
Capitalize "left" and "right" as "LEFT" and "RIGHT":	<input type="checkbox"/>	[False]
Capitalize "living will" as "LIVING WILL":	<input type="checkbox"/>	[False]
Expand and capitalize "NKDA" as "NO KNOWN DRUG ALLERGIES":	<input type="checkbox"/>	[False]
Spell "PO2" and "PCO2" as "pO2" and "pCO2":	<input checked="" type="checkbox"/>	[True]
Capitalize "RBC/hpf" and "WBC/hpf" as "RBC/HPF" and "WBC/HPF":	<input type="checkbox"/>	[False]
Capitalize "trauma alert" as "TRAUMA ALERT":	<input checked="" type="checkbox"/>	[True]
The words "room" and "emergency room" will be capitalized before a number:	<input checked="" type="checkbox"/>	[True]

[Check All](#) | [Uncheck All](#) | [Restore Defaults](#)

Bracketed text in green indicates default values.
Text in bold indicates settings whose value has changed from default.

Person Name Capitalization

<input checked="" type="checkbox"/>	Selection	Results
* = Default selection		
<input type="checkbox"/>	Standard*	Jason T. Myers
<input type="checkbox"/>	All caps	JASON T. MYERS

Department Name Capitalization

<input checked="" type="checkbox"/>	Selection	Results
<input type="checkbox"/>	Standard*	Nuclear Medicine Department
<input type="checkbox"/>	All caps	NUCLEAR MEDICINE DEPARTMENT

*** = Default selection**

Drug Name Capitalization

<input checked="" type="checkbox"/>	Selection	Results
<input type="checkbox"/>	Initial caps	Aspirin Buspar
<input type="checkbox"/>	Standard*	Aspirin BuSpar
<input type="checkbox"/>	All caps	ASPIRIN BUSPAR

*** = Default selection**

Capitalize the Word After a New Line

<input checked="" type="checkbox"/>	Selection	Results
<input type="checkbox"/>	Checked*	First word after a line break is automatically capitalized.
<input type="checkbox"/>	Not checked	First word after a line break is not automatically capitalized, unless another capitalization rule applies to it.

*** = Default selection**

Capitalization/Expansion of Special Words and Phrases

<input checked="" type="checkbox"/>	Selection	Results
<input type="checkbox"/>	allergy, allergies	ALLERGY, ALLERGIES
<input type="checkbox"/>	against medical advice*	AGAINST MEDICAL ADVICE
<input type="checkbox"/>	AMA	AGAINST MEDICAL ADVICE
<input type="checkbox"/>	code blue*	CODE BLUE
<input type="checkbox"/>	code R*	CODE R
<input type="checkbox"/>	code status*	CODE STATUS
<input type="checkbox"/>	code S*	CODE S
<input type="checkbox"/>	code T*	CODE T
<input type="checkbox"/>	do not intubate do not resuscitate*	DO NOT INTUBATE DO NOT RESUSCITATE
<input type="checkbox"/>	DNI DNR	DO NOT INTUBATE DO NOT RESUSCITATE
<input type="checkbox"/>	durable power of attorney*	DURABLE POWER OF ATTORNEY
<input type="checkbox"/>	fish	FISH
<input type="checkbox"/>	left and right	LEFT, RIGHT
<input type="checkbox"/>	living will	LIVING WILL
<input type="checkbox"/>	NKDA	NO KNOWN DRUG ALLERGIES
<input type="checkbox"/>	PO2 PCO2*	pO2 pCO2
<input type="checkbox"/>	RBC/hpf WBC/hpf	RBC/HPF WBC/HPF
<input type="checkbox"/>	trauma alert*	TRAUMA ALERT
<input type="checkbox"/>	room; emergency room*	ROOM 25; EMERGENCY ROOM 1

* = Default selection