



Dragon Medical Enterprise Edition
Version 10.10.304.107
Release Notes

Copyright

2012. Dragon Medical Enterprise Edition, Version 10.10.304.107, Release Notes

This material may not include some last-minute technical changes and/or revisions to the software. Changes are periodically made to the information described here. Future versions of this material will incorporate these changes.

No part of this material or software may be reproduced in any form or by any means, including, without limitation, electronic or mechanical, such as photocopying or recording, or by any information storage and retrieval systems, without the express written consent of Nuance Communications, Inc. Specifications are subject to change without notice.

Copyright © 2008-2012 Nuance Communications, Inc. All rights reserved.

Nuance, the Nuance logo, Dragon, NaturallySpeaking, NaturallyMobile, RealSpeak, Nothing But Speech (NBS), Natural Language Technology, Select-and-Say, BestMatch, MouseGrid, and Vocabulary Editor are registered trademarks or trademarks of Nuance Communications, Inc. All other names and trademarks referenced herein are trademarks of Nuance Communications, Inc., or their respective owners.

Release 10.10.304.107

TTP defects fixed in Dragon Medical Enterprise release 10.10.304.107

For more details, please see search for the TTP number in Test Track Pro web at <http://10.1.4.41/ttweb/login.htm>

Note: This document only includes defects fixed since the last released version of Dragon Medical Enterprise, version 10.10.300.073

TTP or Notes Bug Number	Description	Request Type	For more information
TTP 544	Allscripts Unstructured Note -- Non Standard Window	Defect	Fixed this problem: Dragon Medical Version 10.10.300.021 Allscripts Enterprise version 11.2 Both applications are installed locally When dictated into the Unstructured Note area - receive message "Non Standard Window"
MIGRO3594 TTPs 307 and 308	Jumping cursor issues in PS360 for user that use text streaming (ported from N6_IRR)	Defect	See bug report. For N5_OEM, this fix was included in the 10.1 US Med SDK (for MedQ) and 10.1 UK Med SDK (for Big Hand)
TTP 323 SEJAN521	Error appears while trying to say anything in Adobe Reader X	Defect	See bug report
TTP 403/ OSMAY369	Dragon Client Update for Citrix Patch (vddnspatch.exe) should be updated to not validate ICA ini parameters	Defect	See bug report
TTP 531/ JABYR14	GE Centricity PACS-IW 3.7.3.7 image viewer disables transfer by Dragon Medical Enterprise Edition of recognized text into any dictation destination	Defect	
ANOGO417, ANOGO418	Enabled XenApp5 and XenApp6 support for Citrix Legacy	Feature Request	
No bug number	Updated Admin Guide with information about using Remote Desktop 7.x with Windows Server 2008 R2	Feature Request	