



Using Dragon Medical Practice Edition on devices with high DPI displays



Copyright

2013. Dragon Medical Practice Edition

Copyright © 2002-2013 Nuance Communications, Inc. All rights reserved.

Nuance, ScanSoft, the Nuance logo, the Dragon logo, Dragon, DragonBar, NaturallySpeaking, NaturallyMobile, RealSpeak, Nothing But Speech (NBS), Natural Language Technology, Select-and-Say, MouseGrid, and Vocabulary Editor are registered trademarks or trademarks of Nuance Communications, Inc. in the United States or other countries. All other names and trademarks referenced herein are trademarks of Nuance Communications or their respective owners. Designations used by third-party manufacturers and sellers to distinguish their products may be claimed as trademarks by those third-parties.

L-3660

Known issues using Dragon Medical Practice Edition on devices with high DPI displays

The following information applies to devices, such as Windows tables and Ultrabooks, where the display magnification is set to 150%.

- NLP commands for some applications that use absolute button or menu locations may not work. For example, Dragon's built-in integration with Outlook 2013.
- The DragonBar overlaps with or is overlapped by maximized application windows. To resolve this issue, you can set the display magnification to 100%.
- When using Dragon in tray-icon mode, you cannot click on the menu items after right-clicking on the tray icon. You have to use voice commands to restore the DragonBar. This only occurs if a user is uses a finger to touch the screen. If the user uses a high resolution pointing device, such as a keyboard or the pen that ships with the Microsoft Surface Pro tablet, the menu is functional.

Please be advised that Nuance does not currently support using Dragon Medical Practice Edition on devices with high DPI displays. You may do so if this system configuration works in your environment.